



DVS Sciences' Reagents Roadmap: Near and Far

May 2013


Scott D. Tanner

Chief Technology Officer

in which I take credit for the brilliance of the R&D team

MaxPar[®] Panel Kits: Key Benefits

- MaxPar Panel Kits are the first complete reagent solution for CyTOF[®] and CyTOF[®] 2 experiments
- Kit includes all the reagents require to conduct 16-17 marker assay
- Panel Kit Contents
 - MaxPar[®] Metal Conjugated Antibodies (MCA) 16-17 per kit
 - MaxPar[®] Cell Staining Buffer
 - MaxPar[®] Fix and Perm
 - MaxPar[®] Intercalator (Ir)
 - MaxPar[®] Water



MAXPAR[®] Panel Kits

● Human Kits

- PB Phenotyping 17 marker
- T cell Phenotyping 16 marker
- PBV Basic Phenotyping 7 marker

● Mouse Kits

- Sp/LN Phenotyping 16 marker
- Sp/LN Basic Phenotyping 6 marker

● additional panels being developed for subsequent availability

sensitivity



Necessary characteristics for effective antibody tag

- stable, tightly bind, maintain ab avidity, ...

Desirable characteristics

- maximum # of tag atoms, universal construct ...

How we approach those

- monodisperse metal-chelating polymer


Currently offer two (effectively equivalent) polymers

- X8
- DN3

sensitivity

Have never shown the structure of the X8 polymer


- but you may have seen a related structure published in 2007:



Lou et al, Angewandte Chemie 46, 6111 (2007)

- X8: different backbone; $n \sim 22$
- DN3: based on a cystamine-core dendrimer; $n \sim 16$

chelators



sensitivity

New tagging polymer in development


● variants on “DM48M”

| Item | Sample | Copies of chelator | Atoms/Ab | Tags/Ab |
|------|------------------|--------------------|----------|---------|
| 1 | DM48M_CD20-147Sm | 100 | 365.89 | 3.7 |
| 2 | DM48M_CD45-154Sm | 100 | 382.40 | 3.8 |
| 3 | DM48M_CD25-149Sm | 100 | 330.67 | 3.3 |
| 4 | DM48I_CD20-147Sm | 100 | 418.38 | 4.2 |
| 5 | DM48I_CD45-154Sm | 100 | 260.32 | 2.6 |
| 6 | DM48I_CD25-149Sm | 100 | 318.43 | 3.2 |
| 7 | X8_CD20-147Sm | 22 | 95.32 | 4.3 |
| 8 | X8_CD45-154Sm | 22 | 84.60 | 3.8 |
| 9 | X8_CD25-149Sm | 22 | 68.56 | 3.1 |

DM48M

X8

Titration for CD20-tag on Ramos



sensitivity

long term goal

- quantum step in sensitivity using nanocrystal probes
 - e.g., a 10 nm dia La nanocrystal \sim 10,000 La atoms
 - poster #B193: Pengpeng et al
 - can make them
 - can functionalize them
 - need to deal with non-specific background issues

new metals

mid- to longer- time frame

- increased mass range facilitates use of ^{89}Y and ^{209}Bi as additional antibody probes
 - both work with current polymer tags
 - Y is particularly low mass = low sensitivity (for “bright” antigens)
 - Bi loading requires special protocol = to be offered in conjugate ab's only
- ^{115}In loading works – expect to offer in conjugated ab's only
- continue to work on introducing new metals to expand the arsenal (prospects include Ru, Hf, Zr)
- isotopically-enriched cisplatin viability probe
- isotopically-enriched Pd isotopes for bar-coding

universality



linkers

- currently use maleimide
 - works well for most human and mouse ab's

for special cases, working on

- amine linkers (sulfo-SMCC)
- azide linkers

beads - passport

● 4-element “calibration” beads

● Ce, Eu, Ho, Lu

● ^{140}Ce (88.5%), ^{142}Ce (11.1%)

● ^{151}Eu (47.8%) ^{153}Eu (52.2%)

● ^{165}Ho (100%)

● ^{175}Lu (97.4%)

● Ce is unique identifier - easy to distinguish

● Ce and Ho oxides inform on plasma temperature

● together define mass response curve = quantification

● after tuning, confirm operational characteristics (=passport)

● added to sample, act as internal standard for quantification


● added to sample, provide basis for normalization

beads - magnetic

informational purposes

- we are acquiring magnetic beads from various manufacturers
 - will determine elemental composition
 - to advise which may be used without introducing lanthanide (or other) contamination
- we will share results so as to help advise which beads are most useful (and if appropriate advise on special precautions)

the evolution of mass cytometry



Discover More. Imagine More.

CyTOF[®] 2 upgrade path



| Specification | | CyTOF | CyTOF 2 | Benefit | Upgrade Path | |
|---------------------|---------------------------------------|------------------------|------------------------|---|----------------------|-----------|
| Sample | Type | Single cell suspension | Single cell suspension | No change | No change | |
| | Introduction - | manual | Syringe, Single loop | Auto-charge syringe, Dual loop | Improves ease of use | Yes |
| | | automatic | 96 well autosampler | 96 well autosampler | No change | No change |
| | Throughput - | events/sec | 1000 | 1000 | No change | No change |
| | | parameters/sec | 30,000 | 30,000 | No change | No change |
| replicate sample CV | | < 3% | < 3% | No change | No change | |
| Hardware | Channels | up to 91 | up to 120 | Increases possible parameters per cell | Yes | |
| | Instrument Response (counts/pg Tb) | 200,000 | 400,000 | Improves detection of low abundance targets | Yes | |
| | Dynamic Range (orders of magnitude) | 4 | 4.5 | Allows wider range of signal detection | Yes | |
| | Tuning, Calibration and normalization | manual | automated | Improves data reproducibility and ease of use | Yes | |