

Multiplexing: Measure more than cytokines!



AGENDA

- ✓ Module 1: xMAP technology based genomic and proteomic applications
- ✓ Module 2: Lab demo of SNP typing assay

- Module 3: Milliplex MAP phospho-signaling**

- Module 4: Lab demo of a signaling assay
- Module 5: Wrap up, additional questions/review

Ramsey McIntire, PhD

Multiplex & Cytometry Specialist

EMD Millipore

MILLIPIX MAP - Market Leader!


Human

- Cytokine/chemokine
- T lymphocytes
- Cytokine receptors
- Immunoglobulin
- Cardiovascular
- Cancer
- IGFBPs
- Angiogenesis
- Metastasis
- Neurodegenerative
- Neuropeptide
- Adipokine/Adipocyte
- Metabolic hormone
- Liver
- Pituitary
- Thyroid
- Bone
- Skin
- Kidney Toxicity



Signaling

- Y/S/T phosphorylation
- Apoptosis
- RTK
- OxPhos
- Fatty Acid Ox
- Glycolysis



Primate

- Cytokine/chemokine
- Metabolic hormone



Rat

- Cytokine/chemokine
- Immunoglobulin
- Adipokine/Adipocyte
- Metabolic hormone
- Pituitary
- Stress hormone
- Thyroid
- Bone
- Cardiovascular
- Neuropeptide
- Kidney toxicity
- Vascular injury



Mouse

- Cytokine/chemokine
- T lymphocytes
- Immunoglobulin
- Adipokine/adipocyte
- Metabolic hormone
- Pituitary
- Thyroid
- Bone
- Cardiovascular
- Angiogenesis
- Neuropeptide



Canine

- Cytokine/chemokine
- Metabolism
- Pituitary
- Kidney Toxicity



Porcine

- Cytokine/chemokine




The Challenges of Cell Signaling

Cell signaling is part of a complex system of communication that dictates basic activities in the cell and coordinates cellular functions.

Protein targets of interest are often transient and low expression.

Why multiplex?

- Analyze multiple proteins in a pathway simultaneously
- Analyze multiple pathways in one cell/tissue lysate simultaneously
- Save samples and time!



Milliplex vs. Western Blot

10 proteins
40 samples

.....

Number of assays required

Total time to result

Results per assay

Total sample used

Internal controls possible?

xMAP[®] Technology




18 hours

400

50 μ l (1-25 μ g)

YES

Western blot




Days/weeks

~10


~200 μ l (~100-500 μ g)

NO

Western Blot Comparison




Milliplex Intracellular Assay Procedure



Akt/mTOR Panel: 11-plex

Analytes: p70S6K (Thr412)
IRS1 (Ser312)
GSK3 α (Ser21)
GSK3b (Ser9)
Akt (Ser473)
PTEN (Ser380)
IR (Tyr1162/Tyr1163)
IGF1R (Tyr1135/Tyr1136)
RPS6 (Ser235/Ser236)
TSC2 (Ser939)
mTOR (Ser2448)



The MILLIPLEX MAP Human Akt/mTOR Phosphoprotein Magnetic Bead Panel 11-plex, is used to detect changes in phosphorylated p70S6K (Thr412), IRS1 (Ser312), GSK3 α (Ser21), GSK3b (Ser9), Akt (Ser473), PTEN (Ser380), IR Tyr1162/Tyr1163), IGF1R (Tyr1135/Tyr1136), RPS6 (Ser235/Ser236), TSC2 (Ser939), and mTOR (Ser2448) in cell lysates using the Luminex® system. The detection assay is a rapid, convenient alternative to Western Blotting and immunoprecipitation procedures. Each kit has sufficient reagents for one 96-well plate assay.

Akt/mTOR Panel

11-plex Akt/mTOR Panel Analysis of Insulin Treated HepG2 and IGF-1 Treated MCF-7 Cells




Figure 1. Multiplex analysis of HepG2 and MCF-7 cells treated with insulin or IGF-1. HepG2 cells stimulated with 10 $\mu\text{g}/\text{mL}$ of insulin (15 min) or MCF-7 cells stimulated with 50 ng/mL IGF-1 (10 min) were assayed. The cells were lysed in MILLIPLEX MAP Lysis Buffer containing protease inhibitors. 20 μg total protein of each lysate diluted in MILLIPLEX MAP Assay Buffer 2 were analyzed according the Assay protocol (lysate incubation at 4 $^{\circ}\text{C}$ overnight). The Median Fluorescence Intensity (MFI) was measured with the Luminex[®] system. The figures represent the average and standard deviation of three replicate wells.

Akt/mTOR Panel: Western to Panel Comparison

Lane#	1	2	3	4	5	6	7	8	9
<u>IR pY1162/1163</u>									
MILLIPLEX [®] MFI	45	42	63	1283	861	1220	1112	843	1467
<u>Akt pS473</u>									
MILLIPLEX [®] MFI	21	8	177	1101	28	776	1111	23	1513
<u>GSK3β pS9</u>									
MILLIPLEX [®] MFI	191	29	309	474	184	429	402	273	128

1. HeLa untreated
2. HeLa treated with lambda phosphatase
3. HepG2 untreated
4. HepG2 + insulin
5. HepG2 + 0.1 μM wortmannin + insulin *PI3K inh.
6. HepG2 + 0.1 μM rapamycin + insulin
7. HepG2 + 10 μM U0126 + insulin
8. HepG2 + 50 μM LY-294-002 + insulin *PI3K Inh.
9. HepG2 + 1 μM Ro-31-8220 + insulin *Rsk2 inh.



Home

Welcome,

Instructions

Create a new batch, view installed protocols, perform a daily activity, or access system information and reports.



Click to Create a new Batch from a new Protocol



Click to Create a new Batch using the highlighted Protocol below

Installed Protocols

Name	Version	Manufacturer	Date
Milliplex pAkt-mTOR 11-plex	1	Millipore	10/22/2015 3:03 PM
Equine Cytokine Alpha2	1		8/28/2015 1:07 PM
Mouse Cytokine	1	EMD	8/4/2015 12:18 PM
Equine Cytokine Alpha Kit - AH	1	EMD	3/25/2015 4:19 PM
TGF Beta Kit	3		1/28/2015 10:29 AM
Gilger cytokines	4		1/28/2015 10:27 AM
Feline Premixed Kit	3		1/28/2015 10:22 AM
Canine 6-Plex_6stdpts	2		1/28/2015 10:20 AM
Canine 13 plex protocol	2		1/28/2015 10:14 AM
Porcine Premixed Cytokines	1		7/15/2014 1:20 PM

Scroll



View



Daily Activities



System Initialization

Fluidics prep



Shutdown

Instrument shutdown routine



Probe and Heater

Adjust the probe, set up the heater.



Command:
System State: Complete
Friday 10/23/2015 12:24 PM

Complete

Stop Pause



Drive Fluid Level:

Waste Fluid Level:

Delta Cal Temp: +0.5°C

XY Status: RA1, 24.5°C



- Current Batch
- Saved Batches
- LIS Results
- Reports



Now displaying batch "Milliplex pAkt-mTOR 11-plex".

Instructions View current run statistics and analyte progress per well. For more analysis and result options for this batch, select Details or go to Saved Batches when batch acquisition completes.

Results

Statistic: Median Analyte: Select an Analyte Current Well: 1,H2 Single Step:

Well	Sample	Run Sta...	b-tubulin	GSK3b	IGF1R	IRS1	Akt	mTOR	p70S6K
1,A1	Background0	Ok	82	30	42.5	31	127	73	544
1,B1	Background0	Ok	93	35	50	36	147	78	497
1,C1	Unknown1	Ok	1203	101	69	155	154	437.5	180
1,D1	Unknown1	Ok	1201.5	100	68	152	155	437	178



Date	Message	Code
10/23/2015 12:17:15 PM	Rinse started at RD1.	54B
10/23/2015 12:17:20 PM	Rinse completed at RD1.	54D
10/23/2015 12:17:21 PM	Clean started at RC1.	54B

Save Image
 Chg. Vol
 Progress



Home



Samples



Batches



Results



Protocols



Maintenance




Admin

- Current Batch
- Saved Batches
- LIS Results
- Reports

Now displaying batch "Milliplex pAkt-mTOR 11-plex".

Instructions

View current run statistics and analyte progress per well. For more analysis and result options for this batch, select Details or go to Saved Batches when batch acquisition completes.



Date	Message	Code
10/23/2015 12:17:15 PM	Rinse started at RD1.	54B
10/23/2015 12:17:20 PM	Rinse completed at RD1.	54D
10/23/2015 12:17:21 PM	Clean started at RC1.	54B

Save Image

Chg. Vol

Progress

System Status

Command: Clean at RC1

System State: Active

Friday 10/23/2015 12:20 PM

Batch Running █

Stop

Pause

Eject

Drive Fluid Level: ✔

Waste Fluid Level: ✔

Delta Cal Temp: +0.6°C

XY Status: RC1, 24.5°C

Power Off



- Current Batch
- Saved Batches
- LIS Results
- Reports

Now displaying batch "Milliplex pAkt-mTOR 11-plex".

Instructions View current run statistics and analyte progress per well. For more analysis and result options for this batch, select Details or go to Saved Batches when batch acquisition completes.

Results

Statistic: Median Analyte: Select an Analyte

Current Well 1,H2 Single Step

Well	Sample	Run Sta...	b-tubulin	GSK3b	IGF1R	IRS1	Akt	mTOR	p70S6K
1,A1	Background0	Ok	82	30	42.5	31	127	73	544
1,B1	Background0	Ok	93	35	50	36	147	78	497
1,C1	Unknown1	Ok	1203	101	69	155	154	437.5	180
1,D1	Unknown1	Ok	1201.5	100	68	152	155	437	178
1,E1	Unknown2	Ok	5514	933.5	2102	1383	194	1889	1182
1,F1	Unknown2	Ok	5614	844	1934	990	123	1778.5	1060
1,G1	Unknown3	Ok	6082.5	469	1410.5	1053	146	1322.5	936
1,H1	Unknown3	Ok	5888	405.5	1402	1225.5	162	1215	876
1,A2	Unknown4	Ok	6489	281	805	1012.5	155	916	813
1,B2	Unknown4	Ok	6722	281.5	851	1064	146.5	981	863
1,C2	Unknown5	Ok	3682	2989	1969	1828	188	883	1943

Date	Message	Code
10/23/2015 12:17:15 PM	Rinse started at RD1.	54B
10/23/2015 12:17:20 PM	Rinse completed at RD1.	54D
10/23/2015 12:17:21 PM	Clean started at RC1.	54B

Save Image

Chg. Vol

Progress

System Status

Command: Clean at RC1

System State: Active

Friday 10/23/2015 12:20 PM

Batch Running

Stop

Pause

Eject

Drive Fluid Level: ✔

Waste Fluid Level: ✔

Delta Cal Temp: +0.6°C

XY Status: RC1, 24.6°C

Power Off

Tissue-Specific Responses of IGF-1/Insulin and mTOR Signaling in Calorie Restricted Rats


Naveen Sharma¹, Carlos M. Castorena¹, Gregory D. Cartee^{1,2,3*}

¹ Muscle Biology Laboratory, School of Kinesiology, University of Michigan, Ann Arbor, Michigan, United States of America, ² Department of Molecular and Integrative Physiology, University of Michigan, Ann Arbor, Michigan, United States of America, ³ Institute of Gerontology, University of Michigan, Ann Arbor, Michigan, United States of America

Abstract

Moderate calorie restriction (CR) (~60% of ad libitum, AL, intake) has been associated with numerous favorable physiological outcomes in many species, and the insulin/IGF-1 and mTOR signaling pathways have each been proposed as potential mediators for many of CR's bioeffects. However, few studies have assessed the widely held idea that CR induces the down-regulation of the insulin/IGF-1 and/or mTOR pathways in multiple tissues. Accordingly, we analyzed the phosphorylation status of 11 key signaling proteins from the insulin/IGF-1 (IR^{Tyr1162/1163}, IGF-1R^{Tyr1135/1136}, IRS-1^{Ser312}, PTEN^{Ser380}, Akt^{Ser473}, GSK3 α ^{Ser21}, GSK3 β ^{Ser9}) and mTOR (TSC2^{Ser939}, mTOR^{Ser2448}, P70S6K^{Thr412}, RPS6^{Ser235/236}) pathways in 11 diverse tissues [liver, kidney, lung, aorta, two brain regions (cortex and cerebellum), and two slow-twitch and three fast-twitch skeletal muscles] from 9-month-old male AL and CR Fischer 344 x Brown Norway rats. The rats were studied under two conditions: with endogenous insulin levels (i.e., AL > CR) and with insulin infused during a hyperinsulinemic-euglycemic clamp so that plasma insulin concentrations were matched between the two diet groups. The most striking and consistent effect of CR was greater pAkt in 3 of the 5 skeletal muscles of CR vs. AL rats. There were no significant CR effects on the mTOR signaling pathway and no evidence that CR caused a general attenuation of mTOR signaling across the tissues studied. Rather than supporting the premise of a global downregulation of insulin/IGF-1 and/or mTOR signaling in many tissues, the current results revealed clear tissue-specific CR effects for the insulin signaling pathway without CR effects on the mTOR signaling pathway.

Measured effects of Calorie Restriction and Insulin treatment in multiple skeletal muscles, brain, liver, kidney, heart, and lung



Overall finding: Calorie restriction showed tissue-specific effects, but more consistently impacted insulin/IGF rather than mTOR

Integrative Genomic and Proteomic Analyses Identify Targets for *Lkb1*-Deficient Metastatic Lung Tumors

Julian Carretero,^{1,2,13} Takeshi Shimamura,^{1,3,13} Klarisa Rikova,⁴ Autumn L. Jackson,⁵ Matthew D. Wilkerson,⁵ Christa L. Borgman,¹ Matthew S. Buttarazzi,^{1,7} Benjamin A. Sanofsky,^{1,7} Kate L. McNamara,^{1,7} Kathleyn A. Brandstetter,^{1,7} Zandra E. Walton,^{1,7} Ting-Lei Gu,⁴ Jeffrey C. Silva,⁴ Katherine Crosby,⁴ Geoffrey I. Shapiro,^{1,3} Sauveur-Michel Maira,⁸ Hongbin Ji,⁹ Diego H. Castrillon,¹⁰ Carla F. Kim,¹¹ Carlos Garcia-Echeverria,⁸ Nabeel Bardeesy,^{1,2} Norman E. Sharpless,^{5,6} Neil D. Hayes,⁵ William Y. Kim,^{5,6} Jeffrey A. Engelman,^{1,2} and Kwok-Kin Wong^{1,3,7,*}


SUMMARY

In mice, *Lkb1* deletion and activation of *Kras*^{G12D} results in lung tumors with a high penetrance of lymph node and distant metastases. We analyzed these primary and metastatic de novo lung cancers with integrated genomic and proteomic profiles, and have identified gene and phosphoprotein signatures associated with *Lkb1* loss and progression to invasive and metastatic lung tumors. These studies revealed that SRC is activated in *Lkb1*-deficient primary and metastatic lung tumors, and that the combined inhibition of SRC, PI3K, and MEK1/2 resulted in synergistic tumor regression. These studies demonstrate that integrated genomic and proteomic analyses can be used to identify signaling pathways that may be targeted for treatment.

Src Family Kinase Panel

Because Phosphoscan and western blot analyses detect a phosphorylation site (Y416) that is common to all SRC family kinases (SFKs) including SRC, FYN, LYN, LCK, HCK, FGR, and YES, we interrogated the phosphorylation status of SFKs affected by LKB1 loss in NCI-H358 and A549 cell lines using a Luminex bead assay.


These changes in SRC and FAK suggest impaired adhesion and increased cellular motility (Yeaman, 2004).



D. Phosphoprotein profile of NCI-H358 cells analyzed by protein arrays.

NCI-H358 cells stably expressing shRNA to *LKB1* (LKB1D) or a non-specific shRNA (NT) were subjected to a kinase array as described in materials in methods.

Increased TGFbeta secretion: EMT?



F. TGF- β release analyzed by multiplex bead assay. *In vitro* secretion of TGFB1, TGFB2 and TGFB3 (below detection limit, not shown) were carried out in NCI-H358 cells stably expressing shRNA to *LKB1* (LKB1D) or a non-specific shRNA (NT). Conditioned media were collected after 24h of cell culture and incubated with beads coated with specific antibodies for TGFB1, TGFB2 and TGFB3 as described in Methods. Data is graphed as mean of 3 replicates (pg/10⁶ cells \pm SD).

MILLIPLEX MAP Cell Signaling Portfolio



Pathway Panels:

Glycolysis
Oxidative Phosphorylation
Oxidative Stress
Pyruvate Dehydrogenase
Akt/mTOR
Early Apoptosis
Late Apoptosis
DNA Damage
Heat Shock Protein
MAPK/SAPK
Mitogenesis RTK
Multi-Pathway
NFkB
Src Kinase
STAT
Human T Cell Receptor
TGFbeta

Phospho-Total 2 Plexes:

Akt/PKB
CREB
Erk/MAPK 1/2
IRS1
JNK
mTOR
p38
STAT3
(all of above can be plexed)

Individual MapMates:

46 and counting!

Thank you!

Ramsey McIntire: Multiplex & Cytometry

Karen Tamul: Field Application Scientist

Michelle Dennis: All things EMD Millipore!