

**LSRFortessa**

	Detector	Mirror	Filter	Fluorochrome Examples	Detector Name in Diva
<b>488</b> Blue	A	690LP	695/40	PerCP-Cy5.5, PerCP-eF710	B 695/40
	B	600LP	670/14	PerCP	B 670/14
	C	505LP	525/50	FITC, AF488, GFP, CFsE, YFP, BB515	B 525/50
	D	X	488/10	SSC	SSC

<b>640</b> Red	A	750LP	780/60	APC-Cy7, APC-H7, APC-AF750, APC-AF780	R 780/60
	B	685LP	710/50	APC-Cy5.5, AF700, AF680	R 710/50
	C	X	670/30	APC, AF647	R 670/30

<b>405</b> Violet	A	750LP	780/60	BV 786, BV785	V780/60
	B	690LP	710/50	BV711, Qdot705	V710/50
	C	630LP	660/20	BV650, Qdot655	V660/20
	D	595LP	610/20	BV605, PO, Qdot625, Qdot605	V610/20
	E	505LP	525/50	BV510, AmCyan, Qdot525, LlvE/DEAD Aqua	V525/50
	F	X	450/50	DAPI, PB, BV421, AF405, LlvE/DEAD violet	V450/50

<b>561</b> Yellow Green	A	735LP	780/60	PE-Cy7	YG 780/60
	B	685LP	710/50	PE-Cy5.5	YG 710/50
	C	640LP	670/30	PE-Cy5, 7AAD	YG 670/30
	D	600LP	610/20	PE-TXrd, AF568, PE-AF610, mCherry	YG 610/20

E	565LP Block	582/15	PE, PKH26, DsRed	YG 582/15
---	-------------	--------	------------------	-----------